


Some users may experience Ethernet download errors when using i-SPEED Control and an i-SPEED TR/3/FS camera. This is mainly due to a clash between machine specification and image cache size which can be resolved by increasing the buffer size in the software.


Connect the camera to the computer as normal and connect to the camera in the i-SPEED Control software.

Now you should be on the Acquire tab and the camera is in the list on the left.

Right click on the image of that camera to get the context menu to appear and then select “Advanced Configuration”



Then Change the Batch Size – BMP from 25 to 50 (This does also effect the jpeg settings that are used for download). Then select the “Set as Default” button so that the setting will be remembered. Then select OK.



This should now solve the Ethernet download issue.